

NATO biodosimetry exercise 2011	<b>GAMMA-H2AX</b>	Page 1 of 2
Date created: 25/05/2011	Protocol for blood sampling, irradiation, storage/transport and processing for gamma-H2AX biodosimetry	Version: 2

<b>1. Aim of the method:</b>	Define standard procedures and conditions for obtaining gamma-H2AX-immunostained lymphocyte samples suitable for biological dosimetry based on microscopic foci scoring.
<b>2. Background:</b>	Immunofluorescence staining of DNA double-strand break-associated nuclear foci in blood lymphocytes enables rapid detection of recent radiation exposure in unstimulated blood lymphocytes.
<b>3. Sample requirements:</b>	
3.1 Nature of sample	Aliquots of 2 ml human blood, X-irradiated and diluted 1:1 with RPMI1640 medium + 10% heat-inactivated FCS <sup>1</sup> ; incubated for 2 h and 24 h at 37°C post irradiation.
3.2 Container:	Blood collection tubes with lithium heparin anti-coagulant
3.3 volume	2 ml blood + 2 ml medium per sample
3.4 Transport conditions:	Ship on wet ice (0-4°C) or with a frozen cold pack, temperature logger and dosimeter chip; follow UN Regulation 650 for packaging
3.5 Actions if requirements are not met:	Request new sample; test old sample, discard if no good. Note the potential for continued signal loss in warm samples.
<b>4. Materials:</b>	
4.1. Equipment and supportive materials:	<ul style="list-style-type: none"> <li>• Cytospin with filter paper and slide holder, e.g. Shandon;</li> <li>• Benchtop centrifuge and 15 ml tubes;</li> <li>• Microcentrifuge and 1.5 ml tubes;</li> <li>• Microscopy slides with surface coating to enhance cell adhesion, e.g. SuperFrostPlus;</li> <li>• Fluorescence/phase contrast microscope with at least x40 lens;</li> <li>• Moist chamber (e.g. lidded plastic box with wet paper);</li> <li>• Coplin jars; micro-pipettors and tips; tissues, cover slips, e.g. 24x60 mm<sup>2</sup>; Parafilm stripes.</li> </ul>
4.2. Reagents:	<ul style="list-style-type: none"> <li>• Phosphate buffered saline (PBS), without Mg &amp; Ca;</li> <li>• Deionized H<sub>2</sub>O;</li> <li>• 0.25% Triton-X100, 0.1% glycine in PBS;</li> <li>• 2% formaldehyde (FA), acid free, in PBS;</li> <li>• Blocking solution (BS; 1xPBS, 1% BSA, 0.1% Tween20);</li> <li>• Anti <math>\gamma</math>-H2AX antibody (e.g. Millipore 05-636, 1:500 in BS);</li> <li>• Fluorophore-conjugated 2<sup>nd</sup> antibody (e.g. AlexaFluor488 goat anti-mouse, Invitrogen A-11029, 1:500 in BS);</li> <li>• Mounting medium with antifade and fluorescent DNA counterstain: e.g. Vectashield with DAPI;</li> <li>• Ficoll-based density gradient cell separation medium for human lymphocyte isolation, e.g. Histopaque 1.077;</li> </ul>
<b>5. Sample processing:</b>	<i>All incubations are at ambient temperature unless stated otherwise</i>
5.1 Isolation of lymphocytes	(1) Layer 4 ml diluted blood carefully <sup>2</sup> onto 4 ml Histopaque in 15 ml tube

NATO biodosimetry exercise 2011	<b>GAMMA-H2AX</b>	Page 2 of 2
Date created: 25/05/2011	Protocol for blood sampling, irradiation, storage/transport and processing for gamma-H2AX biodosimetry	Version: 2

5.2 Cytospin	(2) Spin at 1200 g for 5 min without brake <sup>4</sup> and transfer buffy coat cell layer into a 15 ml tube.
5.2A [Alternative: manual spotting]	(3) Add 10ml PBS and swirl briefly; spin 5 min at 400 g; discard supernatant; repeat step (3) once.
5.3 Fixation and permeabilisation	(4) Resuspend pellet in 2 ml PBS. Adjust cell concentration to cytopsin 60-100 000 cells onto a slide for 1 min at 500 RPM. <sup>5</sup>
5.4 Antibody staining	(4A) [Alternative: Resuspend pellet in 0.1 ml PBS and spread 25 $\mu$ l onto a slide. Allow cells to adhere for ~15 min.] <sup>6</sup>
5.5 Mounting	(5) Inspect slides using a phase contrast microscope for suitable density, even distribution and integrity of cells. <sup>7</sup>
5.6 Notes	(6) Fix cells with 2% FA for 5 min. [Optional: Store slides overnight in 0.5% FA/PBS at 4°C]
	(7) Extract cells in 0.25%Triton-X100, 0.1% glycine/PBS for 5 min <sup>8</sup>
	(8) Drain and incubate slides in BS for 10min.
	(9) Remove excess liquid, apply 90 $\mu$ l $\gamma$ -H2AX Ab solution, cover with Parafilm and incubate in moist chamber for 45 min. <sup>9</sup>
	(10) Remove Parafilm with forceps and wash 3 x 2 min in BS.
	(11) Remove excess liquid, apply 90 $\mu$ l 2 <sup>nd</sup> Ab solution, cover with Parafilm and incubate in moist chamber in the dark for 30 min. <sup>9</sup>
	(12) Remove Parafilm and wash 3 x 2 min in PBS in the dark.
	(13) Drain excess liquid, add 18 $\mu$ l mounting medium with DNA counterstain <sup>10</sup> and apply a cover slip.
	(14) Store slides at 4°C until analysis. <sup>11</sup>
	(15) Score foci per cell in random fields of view, using at least an x40 lens. Record results for 20, 30 and 50 scored cells.
	<sup>1</sup> This dilution step reduces haemolysis and enhances lymphocyte recovery. Heat inactivated FCS prevents lymphocytes reacting to 'complement'.
	<sup>2</sup> Tilt tube and let diluted blood run slowly along the tube wall to avoid any mixing.
	<sup>3</sup> Specific isolation of the T-cell subset may reduce uncertainties in the dose estimate, caused by heterogeneous foci formation in different lymphocyte subsets.
	<sup>4</sup> Active braking at the end of the centrifugation spreads out the buffy coat layer.
	<sup>5</sup> Exact cell numbers, volumes and centrifugation time/RPM depend on the specifications of the cytopsin equipment. Keep remaining cells on ice or at 4°C as backup, in case something goes wrong.
	<sup>6</sup> 25 $\mu$ l suffice to cover a ~100mm <sup>2</sup> area on the slide. Smaller volumes should be used in combination with smaller slide areas, e.g. on multi-well slides. Avoid complete drying of the sample. Keep remaining cells on ice or at 4°C as a backup, in case something goes wrong.
	<sup>7</sup> Especially make sure cells form a monolayer. At high densities cells can form multiple layers on top of each other which complicates foci scoring, especially using unsupervised automated approaches.
	<sup>8</sup> Glycine captures unbound FA, thus eliminating the need for additional washing steps.
	<sup>9</sup> Avoid drying of the cell area, as that would result in non-specific staining.
	<sup>10</sup> DNA counterstain can alternatively be added to the 2 <sup>nd</sup> Ab solution to reduce non-specific background fluorescence.
	<sup>11</sup> For long-term storage, slides should be sealed with nail varnish and stored at 4°C for weeks/months or at -20°C for longer periods. Repeated freezing and thawing should be avoided.