

Expanded View Figures




Figure EV1. The poly-GP immunoassay is reproducible.

A–D Poly-GP sandwich immunoassay with anti-GP antibodies 18H8 and 3F9 was used to analyze the GST-GP₁₅ standard at four concentrations. Background-corrected absolute values, mean, and standard deviation (SD) for $n = 4$ GST-GP₁₅ intra-plate replicates (A), $n = 3$ inter-plate replicates (B), and $n = 3$ day-to-day replicates (C). Mean, SD, and the coefficient of variance (CV) for all conditions are listed in (D).